

報 文

末 長 三 雄


矢 田 公 彦

大 深 協 研

microgram concentration is ... 1.0004 ...

00 (1)

STTA 3 mg と四塩化炭素 200 ml を分液漏斗に取




STTA 3 mg と四塩化炭素 200 ml を分液漏斗に取

Order scale in g/h (mm)

80 100




CCl_4 on the stripping of Hg by CN





(6)

Fig. 5 Effect of co-existing actions



4. 検量線と繰り返し精度

Sample	CN added	CN found	Recovery (%)	CN (ppb)
Sea water A (100ml)	0	0.28	—	2.8
	2.0	2.24	98	
	4.0	4.08	94	
Sea water B (100ml)	0	0.25	—	2.5
	2.0	2.18	97	
	4.0	3.98	92	